## Hot Melt Extrusion Process Optimization and Formulation Development of **Amorphous Solid Dispersions for a Poorly Soluble Calcium Channel Blocker**

#### PURPOSE

- 1. Hot melt extrusion (HME) process technology provides the pharmaceutical industry with a means to manufacture amorphous solid dispersions (ASDs) with advantages, such as solvent-free processing, reduced cost of goods (COGs) and continuous manufacturing, making it an attractive option for thermally stable active pharmaceutical ingredients (APIs).
- 2. A 25:75 Z-160:HPMCAS-MMP ASD formulation was successfully processed using a bench-scale mini conical twin-screw (HAAKE Mini-CTW) extruder (Lee et. al., Poster M1430-04-24, AAPS 2019). This goal of the current study was to scaleup the HME process for the lead ASD using a Pharma 11 twin-screw extruder.

### **OBJECTIVE**

1. To perform initial process range finding studies on a Pharma 11 hot melt extruder followed by a **design** of experiment (DOE) approach, using a series of HME processing parameters for process optimization.

### **METHODS**

- 1. A design of experiment (DOE) approach using a series of HME processing parameters (high/low energy mixing input, temperature of heating zones, and screw rate) was applied for the lead ASD process optimization with minimal chemical degradation.
- 2. All of the ASD extrudates were milled and passed through No. 60 sieve to control the ASD particle size range. ASD was characterized by using X-ray powder diffraction, modulated differential scanning calorimetry, and particle size distribution measurement.
- 3. In vitro dissolution performance of Z-160 ASD was evaluated by a two stage non-sink dissolution test in biorelevant dissolution media (0.1 N HCI and fasted-state simulated intestinal fluid), simulating pH and bile salt concentrations for both gastric and intestinal exposure.

F	R	E		SI
	T	āb	le	1.
				Fo
	i			Баю Lo
	I	S	Scr	rew
		Ba	rre	el Te
	ł			
	i			
	I			
	ł		Fe	ed
	Ì	)	Sc	rew
				Τοι
		Ex	ktr	uda
	1		Т	ab
			N	/M
			е	ne
	2	2.	F	Ϊgι
			u	p ·
			а	nd
			Z	on
			С	na
	3	3.	F	igu
				60 200
			to	5au 7 C
		1	Δ	
		Γ.	d	las
			9 W	/ith
			to	o th
I		1000	۰ ۲	
				2.0
		800		J
		000		
	•	600		
	gA/mL	000	1	
	-160] (µ			
	Z	400	1	
		200	1	۶
		0	0	

#### Hong-Guann Lee, Jonathan Yasosky, Pranjal Taskar, Brian Greco, Justin Hughey, Jamie James, Sanjay Konagurthu

hongguann.lee@thermofisher.com | sanjay.konagurthu@thermofisher.com **Thermo Fisher Scientific** 

### ULTS

HME process summary of 25:75 Z-160:HPMCAS-MMP ASD.

rmulation	25:75 Z-160:HPMCAS-MMP ASD							
tch Size (g)	~25-30g							
ot Number	22-1	22-2	23-2	23-1	24-1	24-2	24-3	24-4
Configuration	1 (30°, 60	(30°, 60°, and 90° kneading element)			2 (30° kneading element only)			
emperature (°C)								
Zone 2	50		50		50		50	
Zone 3	125		125		125		125	
Zone 4	160		170		150		140	
Zone 5	160		170		150		140	
Zone 6	160		160		150		140	
Zone 7	160		160		150		140	
Zone 8	160		160		150		140	
Die	160		160		150		140	
Rate (Kg/hr)		0.25						
v speed (rpm)	250	100	250	100	250	100	250	100
rque (Nm)	~4.0	~4.8	~3.5	~7.0	~5.0	~5.6	~6.2	~7.2
ate Appearance	Transparent yellowish							

ole 1 depicts the HME process summary for 25:75 Z-160:HPMCAS-IP ASDs, showing the HME processing parameters with different ergy mixing inputs, temperature of heating zones, and screw speeds.

ure 1 shows two different screw configurations on a Pharma 11 scaletwin-screw extruder with high energy mixing input (configuration 1) low energy mixing input (configuration 2), demonstrating that mixing nes with different kneading elements would affect physical aracteristics of the scale-up Z-160 ASD.

ure 2 shows overlays of XRPD diffractograms for the 25:75 Z-HPMCAS-MMP ASDs, manufactured by Pharma 11 (Table 1). All d ASDs are amorphous, showing a typical amorphous halo compared crystalline Z-160.

scale-up ASDs of 25:75 Z-160:HPMCAS-MMP HME had a single ss transition temperature ( $T_{a}$ ) in Table 2, indicating good homogeneity hout phase separation. The single  $T_{a}$  of each scale-up ASD is similar the  $T_{a}$  (79°C) of feasibility ASD manufactured from Mini-CTW.



Figure 3. Non-sink dissolution profiles of 25:75 Z-160:HPMCAS-MMP scale-up ASDs.

#### **ADVANCING PHARMACEUTICAL SCIENCES, CAREERS, AND COMMUNITY**





Figure 2. XRPD diffractograms of 25:75 Z-160:HPMCAS-MMP scale-up ASDs. Table 2. MDSC results of 25:75 7-160:HPMCAS-MMP scale-up ASDs

able El med de l'estats of Estats de l'estats de l'e					
<b>Processing Condition</b>	Glass Transition Temperature (°C)				
Low energy mixing/250 rpm at 150°C	77.6				
Low energy mixing/100 rpm at 150°C	79.2				
Low energy mixing/250 rpm at 140°C	80.1				
Low energy mixing/100 rpm at 140°C	80.0				
High energy mixing/250 rpm at 160°C	80.0				
High energy mixing/100 rpm at 160°C	80.5				
High energy mixing/100 rpm at 170°C	80.6				
High energy mixing/250 rpm at 170°C	76.6				

5. Figure 3 shows non-sink dissolution profiles of Pharma 11 ASDs, HAAKE mini-CTW ASDs and crystalline Z-160. ASDs processed by high energy mixing input (screw configuration 1) exhibited lower Z-160 supersaturation and sustainment compared with the feasibility ASD processed from mini-CTW. Increasing processing temperature using screw configuration 1 showed low Z-160 dissolution performance, probably due to polymer degradation.

# Pharm Sci 360

### RESULTS

Table 3. Particle size distribution of 25:75 Z-160:HPMCAS-MMP Pharma 11 ASDs using screw configuration 2.

Processing Condition	Dv10 (μm)	Dv50 (μm)	Dv90 (μm)
Low energy mixing/250 rpm at 150°C	158	250	389
Low energy mixing/100 rpm at 150°C	119	207	336
Low energy mixing/250 rpm at 140°C	125	215	364
Low energy mixing/100 rpm at 140°C	89	168	286

- 6. With the exception of ASD processed at 250 rpm/150°C, most ASDs processed by low energy mixing input (screw configuration) 2) exhibited similar Z-160 supersaturation and sustainment compared with the feasibility ASD processed from HAAKE mini-CTW in Figure 3. Screw configuration 2/250 rpm/150°C ASD exhibited significantly lower dissolution performance compared to the HAAKE mini-CTW ASD, probably due to high screw shear stress at 250 rpm.
- Decreasing the processing temperature using screw configuration 2, improved the scale-up ASD dissolution performance, which matches feasibility ASD dissolution performance, indicating that low energy mixing input at 140°C (Z-160  $T_m \sim 127$ °C) is sufficient to manufacture the 25:75 Z-160:HPMCAS-MMP ASD.
- 8. As shown in Figure 3, three ASDs processed by screw configuration 2 exhibited faster dissolution rate than the ASD prepared on the HAAKE mini-CTW at 20 and 40 minutes after gastric transfer in FaSSIF medium. Table 3 shows the particle size distribution difference of the Pharma 11 milled HMEs, and demonstrates that particle size changed Z-160 initial dissolution rate but it did not have any impact on the Z-160 supersaturation and sustainment.

## CONCLUSIONS

- 25:75 Z-160:HPMCAS-MMP ASD was successfully scaled up to a Pharma 11 twin screw extruder with optimized processing parameters from a HAAKE mini-CTW.
- 2. Processing with low energy input and extrusion temperature facilitated the scale-up of 25:75 Z-160:HPMCAS-MMP ASD.
- 3. Particle size differences post milling had an impact on the initial dissolution rate of Z-160 ASD, indicating that ASD particle size may be critical to dissolution performance.

