Summary of capabilities

Groningen, Netherlands

Facility facts:

Workforce: 200

Regulatory approval: EMA&FDA

Contact info: Zuiderweg 72/2

Groningen, 9744 AP The Netherlands

Tel: +31 (0)50 5222 222

Fax: +31 (0)50 5222 333

Unique offering:

Mammalian cell culture manufacturing. Scale up and cGMP manufacturing of recombinant proteins and monoclonal antibodies.

Specialized capabilities:

- cGMP manufacturing for early and late clinical and commercial supply, including process validation and PPQ
- Over 6 years (>70 batches produced) commercial manufacturing experience
- Wide range of Single-Use technologies
- · Modular set-up to match process needs
- Experience in fed batch, perfusion, n-1 perfusion, and XO® processes
- Experience with all commercial cell lines

Development offerings:

| Analytical development equipment | Assay examples |
|--|---|
| | Size exclusion chromatography |
| | Analytical proteinA chromatography |
| HPLC/UPLC based assay using both UV and fluorescent detectors | Ion exchange chromatography (CEX, AEX) |
| | Glycan analysis |
| | Peptide-mapping (UV) |
| | Product specific binding activity |
| 96-384 well ELISA plate reader with UV/Vis and fluorescent detection | Product identity |
| | Impurities: HCP, ProteinA , Long IGF, total protein |
| Spectrophotometer (UV/Vis) Solo VPE (UV/Vis) | Protein concentration |
| Capillary electrophoresis | R/NR CE-SOS |
| Capillary isoelectric focusing (Ice3) | clEF for identity and purity |
| LAL reader (96 well) | Endotoxins |
| TOG reader | Total organic carbon |
| Bioburden filtration setup | Bioburden |
| FTIR | Raw material ID testing |

For detailed capabilities and capacity information please contact your Thermo Fisher Scientific representative.



Groningen key equipment list:

| Proc | cess | Equipment | |
|--------------------|------------|--|--|
| cGMP manufacturing | Upstream | Seed train | Shake flask Rocker-Mixed bioreactor and up to 250 L Single-Use stirred tank bioreactors |
| | | Fed-Batch | 50 L, 1000 Land 2000 L Single-Use stirred tank bioreactors |
| | | Perfusion | 250L or 500L Single-Use bioreactors (ATF and TFF-technology) |
| | | Clarification | Pall Stax & Merck Millistak+ POD 3M Zeta plus encapsulated system |
| | Downstream | Column chromatography | Prepacked columns Single-Use bioprocess systems (3 Uh -510 L/h) Single-Use bioprocess systems (45 L/h - 3500 L/h) |
| | | TFF filtration | Single-Use automatic TFF systems with 1.14 - 20 m2 membrane area capacity Single-Use TangenXTM SIUSTM holder for 0.5 m2 - 5 m2 vmembrane area capacity |
| | | Bulk drug substance final filtration & filling | Disposable sterilizing grade filtration and filling assemblies |
| | | Virus-Filtration (Nano- Filtration) | Planova Filtration from 1 - 4 m2 Filters Virus solve filtration Single-Use membranes up to 12.5m2 Virosart filtration BioEx filtration |

For detailed capabilities and capacity information please contact your Thermo Fisher Scientific representative.



From molecule to medicine: An integrated partner for every step in your drug development journey.

From the earliest stages of development to final distribution to your patient, together, we can make an impact. With more than 65 facilities across 25 countries, Thermo Fisher Scientific pharma services provides end-to-end capabilities across all phases of development, including API, biologics, viral vectors, cGMP plasmids, formulation, clinical trials solutions, logistics services and commercial manufacturing and packaging. We give pharma and biotech companies of all sizes instant access to technical experts across the Americas, Europe, Asia and Australia. Through our Quick to Care™ program, we provide drug development and clinical services that are tailored to fit your unique drug development journey.



Discover the power of our global network.

